



# Tri Auto ZX2




# **Apical Control**

Observe the File Tip Position with the Apex Locator Function







## **Automatic Stop Safety Function**

The Tri Auto ZX2 motor is linked with the canal measurement function. This feature makes endodontic treatment safer and more efficient. The position of the file tip inside the canal is shown in the display, and the motor stops or reverses rotation as soon as the file reaches the end of the working length to prevent over instrumentation. Also, the auto start and stop function automatically starts the motor when the file is inserted into the canal and stops it when it is removed.


**OAS (Optimum Apical Stop) Function**When the file reaches the end of the

working length, it reverses slightly. This reduces the possibility of file jamming and separation.

### **Precise Manual Canal Measurement**

Canals can also be measured manually. A file holder is provided. The ratio of the impedances of two different frequencies is calculated so that the location of the apex is determined with great accuracy regardless of whether the canal is wet or dry. The file tip location is clearly shown in an easy to read meter display. The combination of tactile feedback and the meter reading makes it easy to determine the location of the file tip. A beeping sound also acts as a guide to its location.










# Apical Patency by Rotary Instrumentation

Now Possible with Rotary Files as Small as #10





A Motor for Apical Patency, Glide Path, and Shaping Tri Auto ZX2 is a compact, cordless motor with a built-in apex locator. It can be used for obtaining apical patency and creating a glide path using rotary instrumentation. Rotary files even as small as #10 may be used for these procedures. Responsive, the handpiece allows subtle and delicate manipulation of the file, corresponding to the movements of an experienced dentist's hand. Tri Auto ZX2 also offers a safe and efficient method for canal shaping.











Clinical images are provided by Dr. Yoshi Terauchi / Dr. Toshihiro Ushikubo / Dr. Tai Gega



Apical Patency and Glide Path Using a Motor **OGP Function** (Optimum Glide Path)


Using #20 or smaller files, the motor can be used to achieve apical patency and create a glide path. It reproduces the subtle and delicate finger movements of an experienced dentist. Even curved or constricted canals can be treated faster and more safely. Both NiTi files sized #20 or smaller, or stainless steel files sized #15 or smaller may be used.




Safer and Efficient Canal Shaping
OTR Function (Optimum Torque Release)

Depending on the load on the file, the motor alternates between forward and reverse rotation with great sensitivity to prevent jamming and file breakage. It is also superior in its ability to follow the shape of even curved canals to reduce ledges and over instrumentation. Tri Auto ZX2 is compatible with standard NiTi files - no need to purchase specialty files.

## Canal Shaping Procedure



### File tip position is monitored in real time.



# File Tip Position Monitored in Real Time

File location and canal length are monitored in real time. Even if the working length is altered by shaping, the display accurately shows the file position with respect to the working length.

# Specifications

Name: Tri Auto ZX2
Model: TR-ZX2

Manufacturer: J. MORITA MFG. CORP.

**Functions**: Root Canal Measurement

OGP (Optimum Glide Path)
OTR (Optimum Torque Release)
OAS (Optimum Apical Stop)

Auto Start / Stop Auto Torque Reverse Auto Apical Reverse / Stop Auto Apical Slow Down Auto Apical Torque Down Auto Torque Slow Down

**Options**: Handpiece Holder, External File

Electrode (with cap), Long File

Holder

Handpiece

Free Running Operation

Speed:  $100 \pm 20 - 1000 \pm 100 \text{ rpm}$ 

Torque: 4 Ncm or more

Battery: Lithium-ion battery (DC 3.7 V)
Weight: Approx. 140 g (including battery

and contra angle)

**Battery Charger** 

Rated Input Voltage: DC 5 V Rated Input Current: 2.4 A

Weight: Approx. 280 g (AC adapter not

included)


**AC Adapter** 

Rated Input Voltage: AC100 - 240 V Rated Input Frequency: 47 - 63 Hz Rated Input Current: 0.4 A

### **Ordering Information**

Tri Auto ZX2 24-5902800

Unit includes: handpiece, contra angle, internal file electrode (pre-installed in the contra angle), guide bur, battery charger, battery, AC adapter, power plugs (1 each of 4 types), probe cord, file holder, contrary electrodes (3), tester, handpiece protective sleeves (30), spray nozzle, maintenance oil.



Industrial Design: f/p design gmbh

## Diagnostic and Imaging Equipment

**Treatment Units** 

Handpieces and Instruments

**Endodontic System** 

Laser Equipment

**Laboratory Devices** 

Educational and Training Systems

Auxiliaries



Development and Manufacturing

#### J. MORITA MFG. CORP.

680 Higashihama Minami-cho, Fushimi-ku, Kyoto 612-8533, Japan T +81. (0)75. 611 2141, F +81. (0)75. 622 4595

## Morita Global Website www.morita.com

Distribution

### J. MORITA CORP.

3-33-18 Tarumi-cho, Suita-shi, Osaka 564-8650, Japan T +81. (0)6. 6380 1521, F +81. (0)6. 6380 0585

#### J. MORITA USA, INC.

9 Mason, Irvine CA 92618, USA T +1. 949. 581 9600, F +1. 949. 581 8811

### J. MORITA EUROPE GMBH

Justus-von-Liebig-Strasse 27a, 63128 Dietzenbach, Germany T +49. (0)6074. 836 0, F +49. (0)6074. 836 299

### MORITA DENTAL ASIA PTE. LTD.

3 Science Park Drive, #01-05 The Franklin Singapore Science Park1, Singapore 118223 T +65. 6779. 4795, F +65. 6777. 2279

### J. MORITA CORP. AUSTRALIA & NEW ZEALAND

Suite 2.05, 247 Coward Street, Mascot NSW 2020, Australia T +61. (0)2. 9667 3555, F +61. (0)2. 9667 3577

### J. MORITA CORP. MIDDLE EAST

4 Tag Al Roasaa, Apartment 902, Saba Pacha 21311 Alexandria, Egypt T +20. (0)3. 58 222 94, F +20. (0)3. 58 222 96

### J. MORITA CORP. INDIA

Filix Office No.908, L.B.S. Marg, Opp. Asian Paints, Bhandup (West), Mumbai 400078, India T  $\pm$ 91-22-2595-3482

### J. MORITA MFG. CORP. INDONESIA

28F, DBS Bank Tower, Jl. Prof. Dr. Satrio Kav. 3-5, Jakarta 12940, Indonesia T +62-21-2988-8332, F + 62-21-2988-8201

### SIAMDENT CO., LTD.

71/10 Mu 5, Thakham, Bangpakong, Chachuengsao 24130, Thailand T +66. 38. 573042, F +66. 38. 573043 www.siamdent.com

Catalog Design: f/p design, Germany

Subject to technical changes and errors.