



Debut

## i-Dixel Web

- Platform independence
- No software installation
- 3D image processing

Thinking ahead. Focused on life.

# Features

## 01. Freedom from platform and simplicity of design



Mac & iPad




Windows

## 02. No software installation

With the latest advancements in web technology, i-Dixel Web gives you the freedom to view your images wherever you want and however you want. Gone are the days of complicated chair-side PC setups, and limited choices of hardware. Mac OS X and even iPads can be used to view and edit data from a Morita X-ray system.


i-Dixel Web runs as a web service on an X-ray server PC included with your Morita X-ray system. It serves as a local and secure web-based dental image processing service that you can access throughout your practice on a wide range of devices.



This product supports Internet Explorer 11 on Windows 7 and Safari 7 on Mac OS X 10.9 and iPad (Retina type) on a client PC.

## 03. 3D image processing

Even though i-Dixel Web is a fully web-based system, the features you've come to expect are not compromised. View Morita's attractive 2D and even 3D images on the device of your choice!



Mac, iPad, Safari, OS X and Mac OS are trademarks of Apple Inc., registered in the U.S. and other countries.

Windows, Windows 7 and Internet Explorer are registered trademarks of Microsoft Corporation in the United States and other countries.

Developed and Manufactured by:

### J. MORITA MFG. CORP.

680 Higashihama Minami-cho, Fushimi-ku,  
Kyoto, 612-8533 Japan

T +81. 75. 611 2141, F +81. 75. 622 4595

[www.morita.com](http://www.morita.com)

### J. MORITA MFG. CORP.

33-18, 3-Chome, Tarumi-cho Suita City,  
Osaka, 564-8650 Japan

T +81. 6. 6380 1521, F +81. 6. 6380 0585

### J. Morita Europe GmbH

Justus-von-Liebig-Str. 27a, 63128 Dietzenbach, Germany  
T +49. 6074. 836 0, F +49. 6074. 836 299

[www.morita.com/europe](http://www.morita.com/europe)

### J. Morita USA, Inc.

9 Mason Irvine, CA 92618, USA

T +1. 949. 581 9600, F +1. 949. 465 1095

[www.morita.com/usa](http://www.morita.com/usa)

### J. Morita Corporation Australia & New Zealand

Suite 2.05, 247 Coward Street, Mascot, NSW 2020,  
Australia

T +61. 2. 9667 3555, F +61. 2. 9667 3577

[www.morita.com/anz](http://www.morita.com/anz)

Subject to technical changes and errors.

JME EN FA 1216 \*2.5