

System Requirements for Mac Computer to Access i-Dixel WEB Server

As of January 21, 2026

(For i-Dixel WEB Version 2.25.1)

1. Introduction


The i-Dixel WEB application can be accessed using a web browser on Windows computers, Mac computers, and tablet devices.

Each accessing device should comply with the system requirements.

This document describes the system requirements for Mac computer.

Please note the following:

- (1) The system requirements are as of January 21, 2026.
- (2) The system requirements are for i-Dixel WEB version 2.25.1.
- * The version of i-Dixel WEB can be viewed by clicking on the top right logo.



- (3) Operation of i-Dixel WEB is not guaranteed for operating systems and web browsers that are not described in this document.
- (4) This document assumes that i-Dixel WEB will operate only within the local intranet.

CAUTION

- Updating to the latest versions of operating systems (Windows, Mac, Tablet) or web browsers may cause unforeseen compatibility issues. Please verify the system requirements before performing any update of the operating system or web browser to ensure trouble-free use and operation.

2. Mac Computer

The requirements for macOS with Safari web browser are as follows:

Table 1. Requirements for macOS and Safari web browser

Browser macOS		Safari			
		15	16	17	18
Monterey	12	✘	✘	✘	✘
Ventura ^{*1}	13		✔	✔	✘
Sonoma	14			✔	✘
Sequoia ^{*2}	15				✘

^{*1}: Support ending soon

^{*2}: Support coming soon

■ : Does not exist

✘ : Does not meet system requirements

✔ : Meets system requirements

CAUTION

- Even if the Operating System and browser version meet these requirements, if the hardware does not meet the recommended configuration, the operability of i-Dixel WEB may be impacted. The suggested minimum hardware configuration is shown below in table 2.

Table 2. Recommended Mac computer configuration

CPU	Intel Core i5 2.5 GHz or newer
RAM	4.0 GB or more