

# System Requirements for iPad to Access i-Dixel WEB Server

As of February 21, 2024  
(For i-Dixel WEB Version 2.21.0)

## 1. Introduction


The i-Dixel WEB application can be accessed using a web browser on Windows computers, Mac computers, and tablet devices.

Each accessing device should comply with the system requirements.

This document describes the system requirements for iPad.

Please note the following:

- (1) The system requirements are as of February 21, 2024.
- (2) The system requirements are for i-Dixel WEB version 2.21.0.
  - \* The version of i-Dixel WEB can be viewed by clicking on the top right logo.



- (3) Operation of i-Dixel WEB is not guaranteed for operating systems and web browsers that are not described in this document.
- (4) This document assumes that i-Dixel WEB will operate only within the local intranet.

## CAUTION

- Updating to the latest versions of operating systems (Windows, Mac, Tablet) or web browsers may cause unforeseen compatibility issues. Please verify the system requirements before performing any update of the operating system or web browser to ensure trouble-free use and operation.

## 2. iPad

Table 1 describes iPad and iPadOS that meet the system requirements.  
The only web browser that has been tested and certified by the factory is Safari.

Table 1. Recommendation of iPad and iPadOS system

iPad <sup>*1</sup> \ iPadOS	Before 15	16 <sup>*2</sup>	17
iPad Air and earlier			
iPad Air 2	×		
iPad (5 <sup>th</sup> generation)	×	✓	
iPad (6 <sup>th</sup> generation)	×	✓	✓
iPad Air (3 <sup>rd</sup> generation)	×	✓	✓
iPad (7 <sup>th</sup> generation)	×	✓	✓
iPad (8 <sup>th</sup> generation)	×	✓	✓
iPad Air (4 <sup>th</sup> generation)	×	✓	✓
iPad (9 <sup>th</sup> generation)	×	✓	✓
iPad Air (5 <sup>th</sup> generation)	×	✓	✓
iPad (10 <sup>th</sup> generation)		✓	✓
iPad Pro 12.9-inch (1 <sup>st</sup> generation)	×	✓	
iPad Pro 9.7-inch	×	✓	
iPad Pro 12.9-inch (2 <sup>nd</sup> generation)	×	✓	✓
iPad Pro 10.5-inch	×	✓	✓
iPad Pro 12.9-inch (3 <sup>rd</sup> generation)	×	✓	✓
iPad Pro 11-inch (1 <sup>st</sup> generation)	×	✓	✓
iPad Pro 12.9-inch (4 <sup>th</sup> generation)	×	✓	✓
iPad Pro 11-inch (2 <sup>nd</sup> generation)	×	✓	✓
iPad Pro 12.9-inch (5 <sup>th</sup> generation)	×	✓	✓
iPad Pro 11-inch (3 <sup>rd</sup> generation)	×	✓	✓
iPad Pro 12.9-inch (6 <sup>th</sup> generation)		✓	✓
iPad Pro 11-inch (4 <sup>th</sup> generation)		✓	✓

<sup>\*1</sup>: Includes both Wi-Fi and Wi-Fi + Cellular models  
How to check iPad model is described in the following Apple web page.

<https://support.apple.com/ja-jp/HT201471>

<sup>\*2</sup>: iPad OS versions 16.1 or later are compatible.  
iPad OS version 16.0 is not supported.

: Does not exist  
× : Does not meet system requirements  
✓ : Meets system requirements

### CAUTION

- Note that the operability of i-Dixel WEB may be significantly impaired due to hardware configuration, installed software, or other factors, even if system requirements are met.